Medical Device Waveform Data **Too Much and Too Fast**

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
Biomedical use case at MGB

Innovation Highlight – Medical Device Waveform Data


Mass General Brigham (ex- Partners Healthcare)

One of the largest Healthcare Provider in US

- # of employee 80,000
- # of Patients 2.5 million

GE Carescape B850

Biomed Use Cases from MGB


- ❖ **Research:** Analysis of blood pressure waveforms to predict hemodynamic parameters such as cardiac output, systemic vascular resistance and aortic compliance
- ❖ **Clinical:** Display of waveforms on a dashboard that also displays alerts based on EMR data


System Overview


Devices

HL7 gateways

Storage

Consumers

Clinical: small data, quick lookback

Quality: medium data, medium lookback

Research: big data, deep lookback

Vendor solutions

High-density storage of high-frequency data, including archive tiers

Mix of apps & tools

A vertical white line with small white circles at the top and bottom, positioned on the left side of the slide.

Technical Challenges

Innovation Highlight – Medical Device Waveform Data

What Are the Actual Messages

```
MSH|^~\&|SMARTLINX_11.4_PHSWEB1511|MGH|||20231116154334.892-  
0500||ORU^R01^ORU_R01|3b16154334901771|P|2.6|||AL|NE||8859/1|||IHE_PCD_001^IHE  
PCD^1.3.6.1.4.1.19376.1.6.1.1.1^ISO|  
PID|||4063657||XXXXXXXXXXXXXXXXX|  
OBR|1|149^device ID002F7AA405E^EUI-64|149^CAPSULE^AC3BCD02F7AA405E^EUI-  
64|69121^MDC_OBS_WAVE_CTS^MDC|||20231116154328.194-0500|20231116154333.194-  
0500||MGHMGHE04I_E0406A|||MGHMGHE04I_E0406A|||DatexA_5.2.19.12_CARESCAPE  
B850_Datex|  
OBX|5|NA|27^^CAPSULE|1.99.2.27|-1^-1^-5^-10^-15^-15^-12^-10^-6^-8^-9^-7^-6^-2^-1^-2 ...  
^-17^-12^-7^-2|uV^uV^UCUM|||F|||20231116154328.194-0500|  
OBX|6|NM|6valuesATTR_TIME_PD_SAMP^MDC|1.99.2.27.1|3.333333|ms^ms^UCUM|||F|||2023  
1116154328  
OBR|2|150^CAPSULE^AC3BCD02F7AA405E^EUI-64|150^CAsampling rate05E^EUI-  
64|69121^MDC_OBS_WAVE_CTS^MDC|||20231116154328.194-0500|20231116154333.194-  
0500||MGHMGHE04I_E0406A|||MGHMGHE04I_E0406A|||DatexA_5.2.19.12_CARESCAPE  
B850_Datex|  
OBX|5|NA|27^^CAPSULE|1.99.3.27|-8^-6^-3^-4^-4^-8^-10^-9^27^25^24^18^ ... 9^-10^-9^-  
7|uV^uV^UCUM|||F|||20231116154328.194-0500|  
OBX|6|NM|67981^MDC_ATTR_TIME_PD_SAMP^MDC|1.99.3.27.1|3.333333|ms^ms^UCUM|||F|||2023  
1116154328.194-0500|
```

start & end time

device ID

values

sampling rate

segment 1

segment 2

Key Challenges


High-velocity, high-volume ingestion

- **Maximize parallelization** to avoid contention points


Variability of device data and metadata

- Support **varying sampling rates** for waveforms, trends and alerts
- **Metadata arrives separately**, often out-of-order with series

Diverse query needs lead to complex storage requirements

- Flexible filter criteria ask for **comprehensive indexing strategy**
- Long retention times for research use case ask for **high-density, tiered storage**


A vertical white line with small white circles at the top and bottom, positioned on the left side of the slide.

Platform Solutions

Innovation Highlight – Medical Device Waveform Data

Partitions, Segments, Buckets & Chunks

SQL Projections

Table: MGB_Test.SegmentReading

Field Name	Datatype
ID	%Library.String
bucketID	%Library.Integer
epoch	%Library.Integer
pos	%Library.Integer
segmentBucket	%Library.String
segmentID	%Library.String
ts	%Library.PosixTime
value	%Library.Integer

Table: MGB_Test.DeviceReading

Field Name	Datatype	C #
ID	%Library.String	
bucketID	%Library.Integer	
deviceID	%Library.String	
epoch	%Library.Integer	
readingID	%Library.Integer	
ts	%Library.PosixTime	
value	%Library.Integer	

Table: MGB_Test.Meta

Field Name	Datatype
ID1	%Library.BigInt
ID	%Library.String
deviceID	%Library.String
epoch	%Library.Integer
patientID	%Library.String
segmentID	%Library.String
facility	%Library.String
unit	%Library.String
room	%Library.String
bed	%Library.String
deviceType	%Library.String
deviceSubId	%Library.String

SQL Queries - Clinical

“Using clinical viewers for viewing patient vital signs along with alarms in close to real-time settings remotely.”

```
SELECT deviceID, deviceType, deviceSubId, ts, value
FROM MGB_Test.SegmentReading sr
JOIN MGB_Test.Meta m ON (sr.segmentID = m.segmentID)
WHERE m.patientID = '0123456789'
      AND sr.ts BETWEEN CURRENT_TIMESTAMP
                AND CAST(CURRENT_TIMESTAMP - 180000000 AS POSIXTIME) -- last 3 minutes
ORDER BY deviceID, ts
```

SQL Queries - Research

“Retrieve all occurrences of arrhythmia alarms for [xyz unit] from the past month, week, or day, along with the corresponding waveforms that occurred immediately before or during each alarm.”

```
SELECT arr.bed, arr.ts, arr.value, ecg.deviceID, ecg.ts, ecg.value
FROM (
    SELECT bed,ts,value
    FROM MGB_Test.SegmentReading sr JOIN MGB_Test.Meta m ON (sr.segmentID = m.segmentID)
    WHERE m.deviceType = 2601 AND sr.value > 3 -- arrhythmia device id and alarm value
    AND m.unit = 'MGHL10' -- specific unit for device
    AND sr.ts >= TO_POSIXTIME('2023-11-01','yyyy-mm-dd')
    AND sr.ts < TO_POSIXTIME('2023-12-01','yyyy-mm-dd') -- November 2023
) arr JOIN (
    SELECT deviceID, bed,ts, value
    FROM MGB_Test.SegmentReading sr JOIN MGB_Test.Meta m ON (sr.segmentID = m.segmentID)
    WHERE m.deviceType = 27 -- ECG device type
) ecg ON ecg.bed = arr.bed -- ECG and arrhythmia devices on same bed
    AND ecg.ts > CAST(arr.ts - 12000000 AS POSIXTIME) -- ECG for 2 min before arr alarm
    AND ecg.ts < CAST(arr.ts + 6000000 AS POSIXTIME) -- and 1 min after
ORDER BY arr.bed, arr.ts, arr.value, ecg.deviceID, ecg.ts
```

A vertical white line with small white circles at the top and bottom, positioned on the left side of the slide.

Test Run in Live Environment

Innovation Highlight – Medical Device Waveform Data

All ICU, OR, and Peri-operative Units at MGH

July 31st 2024 (Wednesday), 12 IRIS ports for all 900-1000 devices

Start	Duration	Total messages	Breakdown	Steady-state
11:12 AM	34.5 min	1,395,399 msgs	96% parameter 4% waveform	710 msgs per second


	MGB Count
Beds	751

Storage Compression

Data set	Message Counts				Stored Size		
	#HL7 messages	By type	#HL7 messages	% of messages	HL7 stored size	IRIS stored size	% HL7 stored size
July 31 st	1395399	Waveform	60339	4%	419	2	0.5%
		Parameter	1335060	96%	8100	999	12%

Time Alignment of Waveform and Parameters

- Location = MGHMGHL06I_L0622A, OBX Type = 133^^CAPSULE
- Location = MGHMGHL06I_L0622A, OBX Type = 27^^CAPSULE
- Location = MGHMGHL06I_L0622A, OBX Type = 3144^^CAPSULE

Result

Ingestion Rate

High-density storage

Efficient clinical & research queries

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