

InterSystems Corporation

IHE-ITI Integration Statement



Vendor	Product Name	Version, Release Date
InterSystems Corporation	HealthShare	2010.1, May 2010
InterSystems HealthShare implements all transactions required in the IHE Technical Framework to support the integration profiles, actors, and options listed below.		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Information Technology Infrastructure Domain (ITI)		
Audit Trail and Node Authentication (ATNA)	Audit Record Repository	
	Secure Application	Encryption
Consistent Time (CT)	Time Client	
Cross Community Access (XCA)	Initiating Gateway	XDS Affinity Domain
	Responding Gateway	
Cross Enterprise Document Sharing (XDS.b)	Document Consumer	
	Document Source	Document Replacement
		Folder Management
	Document Registry	Patient Identity Feed HL7v3
	Document Repository	Document Addendum
		Document Replacement
Document Transformation		
Patient Demographics Query v3 (PDQv3)	Patient Demographics Consumer	
	Patient Demographics Supplier	
Patient Identifier Cross-Reference v3 (PIXv3)	Patient Identity Consumer	Update Notification
	Patient Identity Manager	
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.ihe-j.org
InterSystems Corporation's IHE Statement: www.intersystems.com/healthshare/ihe		

InterSystems Corporation

IHE-PCC Integration Statement



Vendor	Product Name	Version, Release Date
InterSystems Corporation	HealthShare	2010.1, May 2010
InterSystems HealthShare implements all transactions required in the IHE Technical Framework to support the integration profiles, actors, and options listed below.		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Patient Care Coordination Domain (PCC)		
Exchange of Personal Health Record Content (XPHR)	Content Consumer	View
		Document Import
		Section Import
		Discrete Data Import
	Content Creator	
Cross Enterprise Sharing of Medical Summaries (XDS-MS)	Content Consumer	View
		Document Import
		Section Import
		Discrete Data Import
	Content Creator	Referral
		Discharge Summary
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.ihe-j.org
InterSystems Corporation's IHE Statement: www.intersystems.com/healthshare/ihe		

InterSystems Corporation

IHE-RAD Integration Statement



Vendor	Product Name	Version, Release Date
InterSystems Corporation	HealthShare	2010.1, May 2010

InterSystems HealthShare implements all transactions required in the IHE Technical Framework to support the integration profiles, actors, and options listed below.

Integration Profiles Implemented	Actors Implemented	Options Implemented
----------------------------------	--------------------	---------------------

Radiology Domain (RAD)

Cross Enterprise Document Sharing for Images (XDS-I.b)	Document Registry	
--	-------------------	--

Regional IHE Statements

In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.ihe-j.org
--	---	--

InterSystems Corporation's IHE Statement: www.intersystems.com/healthshare/ihe